



anti- EAAT3 antibody

Product Information


Catalog No.:	FNab02606
Size:	100 μ g
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20 $^{\circ}$ C for 24 months (Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	solute carrier family 1(neuronal/epithelial high affinity glutamate transporter, system Xag), member 1
Synonyms:	EAAC1, EAAT3, SLC1A1
Observed MW:	66 kDa
Uniprot ID :	P43005

Application

Specificity:	Human, Mouse, Rat; other species are not tested. Please decide the specificity by homology
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human brain using FNab02606(SLC1A1 antibody) at dilution of 1:50



mouse ovary tissue were subjected to SDS PAGE followed by western blot with FNab02606(SLC1A1 antibody) at dilution of 1:500