



FineTest

Wuhan Fine Biotech Co.,Ltd.

C6-323 Biolake, No. 666 Gaoxin Ave.

East Lake High-tech Development District

Wuhan City, Hubei, China (Zip code 430072)

Tel: 0086-027-87384275 Fax: 0086-027-87800889

HRP conjugated Goat anti Mouse IgM

Product Information

Catalog No.: FNSA-0091

Size: 100 μ 1/500 μ 1

Form: liquid

Conjugation: HRP purified from horseradish root

Purity: The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

Stabilizer: 5mg/ml BSA (IgG-Free, Protease-Free)

Storage: 0.01M Sodium phosphate, 0.25M NaCl, 0.1%proclin300, 50% glycerol, pH7.4.
-20 $^{\circ}$ C for 24 months (Avoid repeated freeze / thaw cycles.).

Application

Tested Application: ELISA, WB

Suggested dilution: 1:2000-1:10000 for western blotting with ECL substrates

1:1000-1:20000 for ELISA and Western blotting with chromogenic substrates

Note: This product is for research use only, not for diagnostic or therapeutic use.