



FineTest

Wuhan Fine Biotech Co.,Ltd.

C6-323 Biolake, No. 666 Gaoxin Ave.

East Lake High-tech Development District

Wuhan City, Hubei, China (Zip code 430072)

Tel: 0086-027-87384275 Fax: 0086-027-87800889

Alexa Fluor 555-Goat Anti-Rabbit IgG(H+L)

Product Information

Catalog No.:	FNSA-0095
Size:	100 μ l
Form:	liquid
Purification:	The antibody was purified by immunoaffinity chromatography
Conjugation:	Alexa Fluor 555
Fluorophore:	Alexa Fluor 555 carboxylic acid
Host:	Goat
Storage:	0.01M Sodium phosphate, 0.25M NaCl, 5mM sodium azide, pH7.6. 4°C, avoid repeated freezing and thawing. Store in dark.

Application

Tested Application:	IF, FC
Recommended dilution:	1:500-1:10000

Note: This product is for research use only, not for diagnostic or therapeutic use.