

anti- ANGPTL2 antibody

Product Information

Catalog No.:	FNab00397
Size:	100 μ g
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20 $^{\circ}$ C for 24 months (Avoid repeated freeze / thaw cycles.)

Background


Induces sprouting in endothelial cells through an autocrine and paracrine action.

Immunogen information


Immunogen:	angiopoietin-like 2
Synonyms:	Angiopoietin like 2, Angiopoietin like protein 2, Angiopoietin related protein 2, ANGPTL2, ARP2, HARP
Calculated MW:	57 kDa
Uniprot ID :	Q9UKU9

Application

Specificity:	Human, Mouse, Rat; other species are not tested. Please decide the specificity by homology
Tested Application:	ELISA, WB, IHC, IP
Recommended dilution:	WB: 1:500-1:5000; IHC: 1:200-1:2000; IP: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human heart using FNab00397(ANGPTL2 antibody) at dilution of 1:50



mouse small intestine tissue were subjected to SDS PAGE followed by western blot with FNab00397(ANGPTL2 antibody) at dilution of 1:500