

anti- Claudin 11 antibody

Product Information


Catalog No.:	FNab01733
Size:	100 μ g
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20 $^{\circ}$ C for 24 months (Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	claudin 11
Synonyms:	claudin 11, CLDN11, OSP, OTM
Calculated MW:	22 kDa
Uniprot ID :	O75508

Application


Specificity:	Human; other species are not tested. Please decide the specificity by homology
Tested Application:	ELISA, WB, IHC, IF, FC
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis using FNab01733(CLDN11 antibody) at dilution of 1:50



Immunofluorescent analysis of HepG2 cells using CLDN11 antibody FNab01733 at 1:100 dilution and PE-labeled goat anti-rabbit IgG(red)



human testis tissue were subjected to SDS PAGE followed by western blot with FNab01733(CLDN11 antibody) at dilution of 1:500