



anti- Inhibin beta B antibody

Product Information


Catalog No.:	FNab04320
Size:	100 μ g
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20 $^{\circ}$ C for 24 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	inhibin, beta B
Synonyms:	Activin beta B chain, INHBB, Inhibin beta B chain, inhibin, beta B
Calculated MW:	45 kDa
Uniprot ID :	P09529

Application

Specificity:	Human, Mouse, Rat; other species are not tested. Please decide the specificity by homology
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemical of paraffin-embedded rat kidney using FNab04320(Inhibin beta B antibody) at dilution of 1:100