

anti- Integrin beta-7 antibody

Product Information


Catalog No.: FNab04355
Size: 100 μ g
Form: liquid
Purification: Immunogen affinity purified
Purity: \geq 95% as determined by SDS-PAGE
Host: Rabbit
IsoType: IgG
Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20 $^{\circ}$ C for 24 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen: integrin, beta 7
Synonyms: Integrin beta 7, integrin, beta 7, ITGB7
Calculated MW: 87 kDa
Uniprot ID : P26010

Application

Specificity: Human, Mouse, Rat; other species are not tested. Please decide the specificity by homology
Tested Application: ELISA, WB, IHC
Recommended dilution: WB: 1:500-1:2000; IHC: 1:20-1:200
Image:



Immunohistochemistry of paraffin-embedded human lymphoma using FNab04355(ITGB7 antibody) at dilution of 1:50



FineTest


Wuhan Fine Biotech Co.,Ltd.

C6-323 Biolake, No. 666 Gaoxin Ave.

East Lake High-tech Development District

Wuhan City, Hubei, China (Zip code 430072)

Tel: 0086-027-87384275 Fax:0086-027-87800889



LO2 cells were subjected to SDS PAGE followed by western blot with FNab04355(ITGB7 antibody) at dilution of 1:1000