

anti- MCTP2 antibody

Product Information

Catalog No.:	FNab05069
Size:	100 μ g
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20 $^{\circ}$ C for 24 months (Avoid repeated freeze / thaw cycles.)

Background


Might play a role in the development of cardiac outflow tract.

Immunogen information


Immunogen:	multiple C2 domains, transmembrane 2
Synonyms:	MCTP2
Calculated MW:	100kd
Uniprot ID :	Q6DN12

Application

Specificity:	Human, Mouse, Rat; other species are not tested. Please decide the specificity by homology
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human brain using FNab05069(MCTP2 antibody) at dilution of 1:100



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab05069(MCTP2 antibody) at dilution of 1:500