

anti- MLK3 antibody

Product Information


Catalog No.:	FNab05215
Size:	100 μ g
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20 $^{\circ}$ C for 24 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	mitogen-activated protein kinase kinase 11
Synonyms:	MAP3K11, Mixed lineage kinase 3, MLK 3, MLK3, PTK1, SPRK
Calculated MW:	60 kDa
Uniprot ID :	Q16584

Application

Specificity:	Human, Mouse, Rat; other species are not tested. Please decide the specificity by homology
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:20-1:200; IF: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer using FNab05215(MLK3 antibody) at dilution of 1:50



FineTest


Wuhan Fine Biotech Co.,Ltd.

C6-323 Biolake, No. 666 Gaoxin Ave.

East Lake High-tech Development District

Wuhan City, Hubei, China (Zip code 430072)

Tel: 0086-027-87384275 Fax:0086-027-87800889



PC-3 cells were subjected to SDS PAGE followed by western blot with FNab05215(MLK3 antibody) at dilution of 1:500