



## anti- PAX3 antibody

### Product Information

Catalog No.:	FNab06169
Size:	100 $\mu$ g
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	$\geq$ 95% as determined by SDS-PAGE
Host:	Rabbit
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20 $^{\circ}$ C for 24 months (Avoid repeated freeze / thaw cycles.)

### Background


Transcription factor that may regulate cell proliferation, migration and apoptosis. Involved in neural development and myogenesis.

### Immunogen information


Immunogen:	paired box 3
Synonyms:	CDHS, HUP2, paired box 3, Paired box protein Pax 3, PAX3, WS1
Calculated MW:	60 kDa
Uniprot ID :	P23760

### Application

Specificity:	Human, Mouse, Rat; other species are not tested. Please decide the specificity by homology
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:500-1:1000
Image:	



IP Result of anti-PAX3 (IP: FNab06169, 4ug; Detection: FNab06169 1:500) with mouse skin tissue lysate 3000ug.



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab06169(PAX3 antibody) at dilution of 1:1000