

anti- PLAGL2 antibody

Product Information

Catalog No.:	FNab06517
Size:	100 μ g
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	$\geq 95\%$ as determined by SDS-PAGE
Host:	Rabbit
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20 $^{\circ}$ C for 24 months (Avoid repeated freeze / thaw cycles.)

Background


Shows weak transcriptional activatory activity.

Immunogen information

Immunogen:	pleiomorphic adenoma gene-like 2
Synonyms:	KIAA0198, PLAGL2, Zinc finger protein PLAGL2
Calculated MW:	55 kDa, 68 kDa
Uniprot ID :	Q9UPG8

Application

Specificity:	Human, Mouse, Rat; other species are not tested. Please decide the specificity by homology
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



Raji cells were subjected to SDS PAGE followed by western blot with FNab06517(PLAGL2 antibody) at dilution of 1:500