

anti- RAB5A-Specific antibody

Product Information

Catalog No.:	FNab07039
Size:	100 μ g
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20 $^{\circ}$ C for 24 months (Avoid repeated freeze / thaw cycles.)

Background


The small GTPases Rab are key regulators of intracellular membrane trafficking, from the formation of transport vesicles to their fusion with membranes. The antibody is specific to RAB5A. It has no cross reaction to RAB5B or RAB5C.

Immunogen information


Immunogen:	RAB5A, member RAS oncogene family
Synonyms:	RAB5, RAB5A, Ras related protein Rab 5A
Calculated MW:	24 kDa
Uniprot ID :	P20339

Application

Specificity:	Human, Mouse; other species are not tested. Please decide the specificity by homology
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:3000; IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human brain tissue slide using FNab07039(RAB5A-Specific Antibody) at dilution of 1:200



human brain tissue were subjected to SDS PAGE followed by western blot with FNab07039(RAB5A-Specific antibody) at dilution of 1:500