

## anti- RBM3 antibody

### Product Information

Catalog No.:	FNab07174
Size:	100 $\mu$ g
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	$\geq$ 95% as determined by SDS-PAGE
Host:	Rabbit
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20 $^{\circ}$ C for 24 months (Avoid repeated freeze / thaw cycles.)

### Background


Cold-inducible mRNA binding protein that enhances global protein synthesis at both physiological and mild hypothermic temperatures. Reduces the relative abundance of microRNAs, when overexpressed. Enhances phosphorylation of translation initiation factors and active polysome formation(By similarity).

### Immunogen information


Immunogen:	RNA binding motif(RNP1, RRM) protein 3
Synonyms:	IS1 RNPL, Putative RNA binding protein 3, RBM3, RNA binding motif protein 3, RNPL
Calculated MW:	17 kDa
Uniprot ID :	P98179

### Application

Specificity:	Human, Mouse, Rat; other species are not tested. Please decide the specificity by homology
Tested Application:	ELISA, WB, IHC, IP
Recommended dilution:	WB: 1:200-1:2000; IP: 1:200-1:1000; IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human testis using FNab07174(RBM3 antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab07174(RBM3 antibody) at dilution of 1:500