

anti- TBRG1 antibody

Product Information

Catalog No.:	FNab08531
Size:	100 μ g
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20 $^{\circ}$ C for 24 months (Avoid repeated freeze / thaw cycles.)

Background

Acts as a growth inhibitor. Can activate p53/TP53, causes G1 arrest and collaborates with CDKN2A to restrict proliferation, but does not require either protein to inhibit DNA synthesis. Redistributes CDKN2A into the nucleoplasm. Involved in maintaining chromosomal stability.

Immunogen information

Immunogen:	transforming growth factor beta regulator 1
Synonyms:	NIAM, TB 5, TBRG1
Calculated MW:	45 kDa
Uniprot ID :	Q3YBR2


Application

Specificity: Human; other species are not tested. Please decide the specificity by homology

Tested Application: ELISA, WB

Recommended dilution: WB: 1:200-1:2000

Image:



293T cells were subjected to SDS PAGE followed by western blot with FNab08531(TBRG1 antibody) at dilution of 1:1000