

Rabbit anti- mouse IgA

Product Information


Catalog No.:	FNab09894
Size:	100 μ g
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20 $^{\circ}$ C for 24 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	mouse IgA
Calculated MW:	60 kDa

Application

Specificity:	Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:1000-1:10000; IHC: 1:2000-1:20000
Image:	



Immunohistochemistry of paraffin-embedded mouse lymph gland tissue slide using FNab09894(mouse IgA antibody) at dilution of 1:5000



FineTest


Wuhan Fine Biotech Co.,Ltd.

C6-323 Biolake, No. 666 Gaoxin Ave.

East Lake High-tech Development District

Wuhan City, Hubei, China (Zip code 430072)

Tel: 0086-027-87384275 Fax:0086-027-87800889



mouse plasma were subjected to SDS PAGE followed by western blot with FNab09894(mouse IgA antibody) at dilution of 1:5000