



FineTest

**Wuhan Fine Biotech Co.,Ltd.**

C6-323 Biolake, No. 666 Gaoxin Ave.

East Lake High-tech Development District

Wuhan City, Hubei, China (Zip code 430072)

Tel: 0086-027-87384275 Fax:0086-027-87800889

## Goat anti- human IgM

### Product Information

Catalog No.: FNab09908  
Size: 100 $\mu$  g  
Form: liquid  
Purification: Immunogen affinity purified  
Purity:  $\geq$ 95% as determined by SDS-PAGE  
Host: Goat  
IsoType: IgG  
Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20 $^{\circ}$ C for 24 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen: human IgM

### Application

Specificity: Human; other species are not tested. Please decide the specificity by homology  
Tested Application: ELISA, WB  
Recommended dilution: WB: 1:1000-1:10000

**Note: This product is for research use only, not for diagnostic or therapeutic use.**