

anti- METTL4 antibody

Product Information


Catalog No.:	FNab09913
Size:	100 μ g
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20 $^{\circ}$ C for 24 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	methyltransferase like 4
Synonyms:	HsT661, methyltransferase like 4, METTL4
Calculated MW:	55 kDa
Uniprot ID :	Q8N3J2

Application

Specificity:	Human, Mouse, Rat; other species are not tested. Please decide the specificity by homology
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500 - 1:2000
Image:	



A549 cells were subjected to SDS PAGE followed by western blot with FNab09913(METTL4 Antibody) at dilution of 1:1000