



Anti- p24 antibody

Product Information

Catalog No.:	FNab09919
Size:	100 μ g
Form:	liquid
Purification:	protein A+G purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Mouse
IsoType:	IgG1
Storage:	0.01MPBS pH 7.4, -20 $^{\circ}$ C for 24 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	Recombinant Capsid protein p24
------------	--------------------------------

Application

Specificity:	HIV-1
Tested Application:	ELISA

FOR RESEARCH USE ONLY!