

## anti- SELPLG antibody

### Product Information


Catalog No.:	FNab09941
Size:	100 $\mu$ g
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	$\geq$ 95% as determined by SDS-PAGE
Host:	Rabbit
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20 $^{\circ}$ C for 24 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	P-selectin glycoprotein ligand 1
Synonyms:	CD162, CLA, PSGL 1, PSGL1, selectin P ligand, SELPLG
Calculated MW:	120 kDa
Uniprot ID :	Q14242

### Application

Specificity:	Human; other species are not tested. Please decide the specificity by homology
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:20 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using FNab09941( SELPLG Antibody) at dilution of 1:100



FineTest


**Wuhan Fine Biotech Co.,Ltd.**

C6-323 Biolake, No. 666 Gaoxin Ave.

East Lake High-tech Development District

Wuhan City, Hubei, China (Zip code 430072)

Tel: 0086-027-87384275 Fax:0086-027-87800889



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab09941(SELPLG Antibody) at dilution of 1:500