

anti- AKT1S1 antibody

Product Information


Catalog No.:	FNab09942
Size:	100μ g
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20℃ for 24 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	Proline-rich AKT1 substrate 1
Synonyms:	AKT1S1, Lobe, PRAS40, AKT1 substrate 1
Calculated MW:	50 kDa
Uniprot ID :	Q96B36

Application

Specificity:	Human, Mouse, Rat; other species are not tested. Please decide the specificity by homology
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:20 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab09942(AKT1S1 Antibody) at dilution of 1:100



FineTest


Wuhan Fine Biotech Co.,Ltd.

C6-323 Biolake, No. 666 Gaoxin Ave.

East Lake High-tech Development District

Wuhan City, Hubei, China (Zip code 430072)

Tel: 0086-027-87384275 Fax:0086-027-87800889



mouse heart were subjected to SDS PAGE followed by western blot with FNab09942(AKT1S1 Antibody) at dilution of 1:1000