



FineTest

Wuhan Fine Biotech Co.,Ltd.

C6-323 Biolake, No. 666 Gaoxin Ave.

East Lake High-tech Development District

Wuhan City, Hubei, China (Zip code 430072)

Tel: 0086-027-87384275 Fax:0086-027-87800889

Rabbit anti- goat IgG(H+L)

Product Information

Catalog No.: FNab09946
Size: 100 μ g
Form: liquid
Purification: Immunogen affinity purified
Purity: \geq 95% as determined by SDS-PAGE
Host: Rabbit
IsoType: IgG
Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20 $^{\circ}$ C for 24 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen: goat IgG
Calculated MW: 155 kDa

Application

Specificity: goat; other species are not tested. Please decide the specificity by homology
Tested Application: ELISA, WB