

anti- ATPB antibody

Product Information


Catalog No.:	FNab09949
Size:	100 μ g
Form:	liquid
Purification:	protein A+G purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Mouse
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20 $^{\circ}$ C for 24 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	ATP synthase, H ⁺ transporting, mitochondrial F1 complex, beta polypeptide
Synonyms:	ATPMB, ATPSB, HEL-S-271
Calculated MW:	56 kDa
Uniprot ID :	P06576

Application

Specificity:	Human, Mouse, Rat; other species are not tested. Please decide the specificity by homology
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:5000 - 1:30000; IHC: 1:200 - 1:1000
Image:	



Immunohistochemistry of paraffin-embedded rat heart using FNab09949(ATPB antibody) at dilution of 1:500



FineTest


Wuhan Fine Biotech Co.,Ltd.

C6-323 Biolake, No. 666 Gaoxin Ave.

East Lake High-tech Development District

Wuhan City, Hubei, China (Zip code 430072)

Tel: 0086-027-87384275 Fax:0086-027-87800889



Hela cells were subjected to SDS PAGE followed by western blot with FNab09949(ATPB antibody) at dilution of 1:10000