

## anti- RAB33B antibody

### Product Information


Catalog No.:	FNab10126
Size:	100 $\mu$ g
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	$\geq$ 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20 $^{\circ}$ C for 24 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	Ras related protein Rab 33B
Synonyms:	RAB33B, Ras related protein Rab 33B
Observed MW:	26 kDa
Uniprot ID :	Q9H082

### Application

Specificity:	Human; other species are not tested. Please decide the specificity by homology
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab10126 (RAB33B antibody) at dilution of 1:1000



**FineTest**

**Wuhan Fine Biotech Co.,Ltd.**

C6-323 Biolake, No. 666 Gaoxin Ave.

East Lake High-tech Development District

Wuhan City, Hubei, China (Zip code 430072)

Tel: 0086-027-87384275 Fax:0086-027-87800889

---